



How to control your device

Wish to turn your NeoRhythm's lights, vibration or Bluetooth off?

Learn how to control your device with tapping, lights and vibration.

1. LIGHTS

- **Program Currently Running:** Green light every second




1 blink per 1 second

- **The device is connected, but no session is running:** Green light every 5 seconds



1 blink per 5 seconds

- **Low Battery:** Red Light every second




1 blink per 1 second

2. VIBRATION

You can turn your vibration settings on in the app.

Connect your device, then tap the settings wheel in the bottom right corner.

Tap the button to turn the vibration ON/OFF.



- **Double vibration:** the session has ended
- **Single vibration:** the device is connected, the session has started

3. TAPPING

The devices are all gesture-controlled. The device can restart the last used session without having to connect to the app first.

- **Double tap:** end the current session
- **Double tap when no session is running:** restart the last used session
- **Single tap:** wakes the device up