

# Light Indicators

## Therapy Modes Overview

NeoSurround offers a variety of therapy programs designed to support your well-being in different areas of life. Each mode is represented by a specific light color, so you can easily recognize which program is active.

● Enhancing Mental Capabilities (Lime)	Supports focus, concentration, and cognitive performance.
● Energy and Vitality (Lime)	Helps you recharge and feel more energized.
● Custom Program (Green)	Allows you to create and save your own program.
● 24 Hours (Green)	Continuous support throughout the day.
● Pain Control (Golden Rod)	Designed to help relieve discomfort and ease tension.
● Sleep 1 (Deep Pink)	Promotes falling asleep more easily.
● Sleep 2 (Deep Pink)	Supports deeper, continuous sleep.
● Lucid Dreaming (Deep Pink)	Supports dream awareness and deep rest.
● Meditation Theta (Light Pink)	Assists in reaching deeper meditative states.
● Meditation Calming and Sync (Light Pink)	Promotes inner balance and calmness.
● Focus Meditation (Light Pink)	Enhances concentration during meditation.
● Mindful Meditation (Light Pink)	Encourages mindfulness and present awareness.
● Quiet Mind (Light Pink)	Helps calm racing thoughts.
● Open Heart (Light Pink)	Promotes emotional openness and balance.
○ Deep Relaxation (Antique White)	Encourages calmness and stress relief.

**Note:** While a program is active, the light on NeoSurround will display a solid color in the program's designated shade.

## System Status Lights

NeoSurround uses different light colors to show its current status:

 Solid blue light	Device is connected, but no program is running.
 Blinking blue light	Device is not connected / ready to connect.
 Solid yellow light	Device is disconnected.
Program color (see manual)	A program is currently running.

## Error Status Lights

If NeoSurround detects an issue, the light will change color to indicate the type of error:

 Red (Fast Blinking)	Error. Please restart the device or contact support.
 Dark Orange (Blinking)	Incorrect power supply is connected. Use the original power adapter.
 Yellow (Slow Blinking)	Warning. The device may not be functioning properly — check connections and restart if needed.